

Technical Information

Rxn-45 Raman spectroscopic probe



Table of Contents

Function and system design 3

Fields of application 3

Laser safety interlock 3

Rxn-45 probe..... 3

Installation 4

Data collection zone: Short 4

Specifications 5

General specifications..... 5

Probe dimensions..... 6

MPE : Ocular exposure..... 6

MPE: Skin exposure 7

Function and system design

Fields of application

The Rxn-45 Raman spectroscopic probe is designed for the needs of bioprocessing pilot and manufacturing sites.

Recommended applications include:

- **Cell culture:** Glucose, lactate, amino acids, cell density, titer, and more
- **Fermentation:** Glucose, glycerol, acetate, methanol, ethanol, biomass, and more

Use of the device for any purpose other than that described, poses a threat to the safety of people and of the entire measuring system and invalidates any warranty.

Laser safety interlock

The Rxn-45 probe, as installed, forms part of the interlock circuit. The interlock circuit is a low-current electrical loop. If the fiber cable is severed, the laser will turn off within 20 milliseconds of the breakage.

NOTICE

Handle probes and cables with care.

Fiber cables should NOT be kinked and should be routed to maintain the minimum bend radius of 6 inches (152.4 mm).

- ▶ Permanent damage may result if cables are not routed appropriately.

The Electro-Optical (EO) fiber cable with its embedded interlock loop must be plugged into the back of the Raman Rxn analyzer for the appropriate channel. The interlock loop is complete when the probe side of the EO fiber cable is plugged into the Rxn-45 probe.

When there is potential for the laser to be energized, the laser interlock indicator light on the probe body is illuminated.

Rxn-45 probe

The Rxn-45 probe with the right angle connection is shown below.




Figure 1: Rxn-45 probe

| # | Name | Description |
|---|---------------------------------|--|
| 1 | Laser interlock indicator light | Illuminated when there is potential for the laser to be energized |
| 2 | Probe tip | Tip of probe for sample interface; 120 mm immersion length |
| 3 | Flange and o-ring | Welded flange and replaceable USP Class VI o-ring to ensure a tight seal with the vessel port/hardware |
| 4 | Captive nut | PG13.5 thread for industry standard sensor housings; welded port connectors available |
| 5 | Fiber optic cable connector | Electro-Optical (EO) fiber connection under spring-loaded fiber connector cap |

Table 1. Rxn-45 probe parts

Installation

During installation, standard eye and skin safety precautions for Class 3B laser products (as per EN-60825/IEC 60825-14) should be observed as described below.

| | |
|------------------|---|
| ⚠ WARNING | <p>Standard precautions for laser products should be observed.</p> <ul style="list-style-type: none"> ▶ Probes should always be capped and/or pointed away from people toward a diffuse target if not installed in a sample chamber. |
| ⚠ CAUTION | <p>The laser input into the probe must not exceed 499 mW.</p> <p>If stray light is allowed to enter an unused probe, it will interfere with data collected from a used probe and may cause calibration failure or measurement errors.</p> <ul style="list-style-type: none"> ▶ Unused probes should ALWAYS be capped to prevent stray light from entering the probe. |
| NOTICE | <p>When installing the probe <i>in situ</i>, the user must provide the strain relief to the fiber optic cable at the probe installation location.</p> |

Data collection zone: Short

All versions of the Rxn-45 probe utilize short data collection zones. The short data collection zone maximizes spectral reproducibility by minimizing the impact of sample opacity, sample color, and transient particulates on the measured Raman spectrum.

Specifications

General specifications

General specifications for the Rxn-45 probe are listed below. Additionally:

Note: Max service pressure ratings do not include the ratings of any fittings or flanges used to mount the probe into the process system. These items need to be independently evaluated and may lower the maximum service pressure of the probe.

| Item | Description | |
|--------------------------------------|---|---|
| Maximum laser power into probe | <499 mW | |
| Relative humidity | Up to 95%, non-condensing | |
| Maximum service pressure (at tip) | ≤ 13.8 Bar (200 psi) | |
| Laser wavelengths | 785 nm and 993 nm | |
| Process connection | PG13.5 thread for industry standard sensor housings; welded port connectors available | |
| Rating | IP-65 | |
| Depth of field | 0.33 mm (FWHM) | |
| Chemical resistance | Limited by materials of construction | |
| Sterilization protocol compatibility | SIP/CIP | |
| Spectral coverage | 785 nm | 150 – 3425 cm ⁻¹ |
| | 993 nm | 200 – 2400 cm ⁻¹ |
| Temperature | Window, at tip | -30°C to +150°C |
| | Probe body | up to +150°C |
| | Temperature ramp | 30°C/min |
| Probe measurements | Immersion length | 120 mm |
| | Diameter | 12 mm |
| | Dimensions | 306 mm x 127 mm x 34 mm |
| Materials of construction | Probe body | 316L stainless steel |
| | Window | Proprietary material, optimized for bioprocesses |
| Wetted, in contact with sample | Adhesive | USP Class VI and ISO993 compatible |
| | Surface finish | Ra 15 with Electropolish to ASME BPE SF4 finish |
| | Fiber optic cable | <ul style="list-style-type: none"> Design: PVC jacketed, proprietary construction Connections: proprietary electro-optic (EO) or FC to EO fiber converter(s) for non-embedded systems |
| Fiber optic cable specifications | Length | 5m (16.4 feet) standard Custom lengths available in 10m increments, starting at 10m |
| | Minimum bend radius | 152.4 mm (6 inches) |
| | Temperature | -40°C to +70°C |
| | Flame resistance | Certified: CSA-C/US AWM I/II, A/B, 80C, 30V, FT1, FT2, VW-1, FT4 Rated: AWM I/II A/B 80C 30V FT4 |

Table 2. General specifications

Probe dimensions

The dimensions for the Rxn-45 probe are shown below.




Figure 2: Rxn-45 probe dimensions

MPE : Ocular exposure

Refer to the tables below from the ANSI Z136.1 standard to calculate the Maximum Permissible Exposure (MPE) for point source ocular exposure to a laser beam.

A correction factor (C_A) may also be required and can be determined below.

| Wavelength λ (nm) | Correction Factor C_A |
|---------------------------|---------------------------|
| 400 to 700 | 1 |
| 700 to 1050 | $10^{0.002(\lambda-700)}$ |
| 1050 to 1400 | 5 |

Table 3. Wavelength dependent correction factor C_A

| Determining Maximum Permissible Exposure (MPE) for Point Source Ocular Exposure to a Laser Beam | | | | |
|---|----------------------------------|-----------------------------------|-----------------------|---|
| Wavelength λ (nm) | Exposure Duration t (s) | MPE Calculation | | MPE where $C_A = 1.4791$ |
| | | (J·cm ⁻²) | (W·cm ⁻²) | |
| 785 and 993 | 10^{-13} to 10^{-11} | $1.5 C_A \times 10^{-8}$ | - | 2.2×10^{-8} (J·cm ⁻²) |
| | 10^{-11} to 10^{-9} | $2.7 C_A t^{0.75}$ | - | Insert time (t) and calculate |
| | 10^{-9} to 18×10^{-6} | $5.0 C_A \times 10^{-7}$ | - | 7.40×10^{-7} (J·cm ⁻²) |
| | 18×10^{-6} to 10 | $1.8 C_A t^{0.75} \times 10^{-3}$ | - | Insert time (t) and calculate |
| | 10 to 3×10^4 | - | $C_A \times 10^{-3}$ | 1.4971×10^{-3} (W·cm ⁻²) |

Table 4. MPE for ocular exposure with 785 nm or 993 nm laser emission

MPE: Skin exposure

Refer to the table below from the ANSI Z136.1 standard to calculate the MPE for skin exposure to a laser beam.

| Maximum Permissible Exposure (MPE) for Skin Exposure to a Laser Beam | | | | |
|--|--------------------------------------|------------------------|-----------------------|--|
| Wavelength λ (nm) | Exposure Duration t (s) | MPE Calculation | | MPE where $C_A = 1.4791$ |
| | | (J·cm ⁻²) | (W·cm ⁻²) | |
| 785 and 993 | 10 ⁻⁹ to 10 ⁻⁷ | 2 $C_A \times 10^{-2}$ | - | 2.9582 × 10 ⁻² (J·cm ⁻²) |
| | 10 ⁻⁷ to 10 | 1.1 $C_A t^{0.25}$ | - | Insert time (t) and calculate |
| | 10 to 3 × 10 ⁴ | - | 0.2 C_A | 2.9582 × 10 ⁻¹ (W·cm ⁻²) |

Table 5. MPE for skin exposure with 785 nm or 993 nm laser emission

www.addresses.endress.com
